



## Fusion® Lateral OA Plus Knee Brace



Sizes: s, m, m+, l, xl, xxl

The Fusion® Lateral OA Plus brace features a streamlined frame with an innovative, slim thumbwheel dial. Designed to provide lateral off-loading for patients with unicompartmental osteoarthritis. The Fusion® Lateral OA Plus incorporates an offset gauge that clearly displays the amount of load applied in half-millimeter increments. It features a flexible polymer, allowing patients to always get an intimate, contoured fit that keeps the brace in place all day long.

[Learn More About Knee Osteoarthritis](#)

### Common Examples of Use

- Unicompartmental OA-lateral
- Chondral defect procedures
- Meniscal procedures

### Product Features

- AirTech® Frame pads - innovative windows and channels in the pad permit cooling air flow while allowing moisture and heat to move away from the skin
- Innovative thumbwheel dial makes precision off-loading adjustments easy - no tools required
- PivotPoint straps provide superior comfort during motion
- Cushioned straps and pads increase comfort
- Sleek frame design provides unhindered mobility

Category: Functional OA. Tag: Fusion Braces.

### Instructions for Use

## Part Numbers

### Part Numbers

<b>Part #</b>	<b>Description</b>
01302	Fusion® Lateral OA Plus Custom
01303	Fusion® Lateral OA Plus Custom, Color/Pattern Option
077XX	Fusion® Lateral OA Plus Prefabricated, Left
079XX	Fusion® Lateral OA Plus Prefabricated, Right
795XX	Fusion® Lateral OA Plus, Refurbish, Left*
797XX	Fusion® Lateral OA Plus, Refurbish, Right*

\* Refurb kit includes 2 stop sets, frame pads with hook and loop fasteners, straps and instructions

## Product Sizing Chart

Fusion Lateral OA Plus Size Chart			
Size #	Size	6" Above Mid Patella	6" Below Mid Patella
XX = 20	S	15" - 18" (38 - 46 cm)	12" - 14" (30 - 36 cm)
XX = 30	M	18" - 19.5" (46 - 50 cm)	14" - 15" (36 - 38 cm)
XX = 35	M+	19.5" - 21" (50 - 53 cm)	15" - 16" (38 - 41 cm)
XX = 40	L	21" - 24" (53 - 61 cm)	16" - 18" (41 - 46 cm)
XX = 50	XL	24" - 27" (61 - 69 cm)	18" - 20" (46 - 51 cm)
XX = 60	XXL	27" - 30" (69 - 76 cm)	20" - 23" (51 - 58 cm)

*Circumference measured 6" (15 cm) above mid patella and 6" (15 cm) below mid patella.  
Not available in XS.  
Brace length starts at 13.5" and varies by size.*