



## Freestyle™ OA Knee Brace



Sizes: s, m, l, xl, xxl

The Freestyle™ OA provides off-loading of the medial or lateral compartment for patients with mild to moderate unicompartmental osteoarthritis. The simple wraparound design allows patients to apply and remove the brace with ease. Comfortable, breathable material provides focused compression around the knee joint. The single-sided low profile hinge encourages patient compliance, while the innovative thumbwheel dial allows precision unloading adjustments without the use of a tool.

[Learn More About Knee Osteoarthritis](#)

### Common Examples of Use

- Mild to moderate unicompartmental osteoarthritis
- Meniscus procedures/repairs
- Osteochondral defect
- Articular cartilage defect/repair

### Product Features

- Easy to apply and remove
- Low profile to provide patient comfort
- Compression provides warmth and decreases swelling
- Innovative thumbwheel dial makes precision off-loading adjustments easy - no tools required

Category: Functional OA.

### Instructions for Use

[Instructions for Use](#)

## Part Numbers

## Part Numbers

<b>Part #</b>	<b>Description</b>
1171X	Freestyle™ OA Medial, Left
1172X	Freestyle™ OA Medial, Right
1176X	Freestyle™ OA Lateral, Left
1177X	Freestyle™ OA Lateral, Right

## Product Sizing Chart



image not found or type unknown