



DUO Knee Brace



Sizes: xs, s, m, m+, l, xl, xxl

A unique design, our DUO (Dynamic Unloading Osteoarthritis) brace is the only dual-upright dynamic OA brace on the market. This brace creates a load across the joint when the knee is extended and reduces the load as the knee flexes. Thus, when the patient is bearing weight, the brace is actively unloading. DUO reduces force when sitting, making it ideal for patients who need to wear a brace for daily activities.

[Learn More About Knee Osteoarthritis](#)

Common Examples of Use

- Unicompartmental OA
- Unicompartmental OA / meniscal instabilities with ligamentous injuries

Product Features

- Low-profile frame
- Dynamic loading hinge
- Incremental arm adjustments
- Up to 16 degrees of correction can be applied
- Unique strapping design minimizes migration

Category: Functional OA.

Instructions for Use

[Instructions for Use](#)

Part Numbers

Part Numbers

Part #	Description
ED312000	DUO Custom Short
ED112000	DUO Custom
ED114YXX--B	DUO, Lateral, XS-XXL
ED112YXX--B	DUO, Medial, XS-XXL
ED113YXX--B	DUO, Lateral Athletic, XS-XXL
ED111YXX--B	DUO, Medial Athletic, XS-XXL
ED314YXX--B	DUO, Lateral, Short, XS-XXL
ED312YXX--B	DUO, Medial, Short, XS-XXL
ED313YXX--B	DUO, Lateral Athletic, Short, XS-XXL
ED311YXX--B	DUO, Medial Athletic, Short, XS-XXL
ED0201XX	Kit, Soft Goods, Left Medial / Right Lateral, XS-XXL
ED0202XX	Kit, Soft Goods, Right Medial / Left Lateral, XS-XXL
ED0211XX	Kit, Soft Goods, Athletic, Left Medial / Right Lateral, XS-XXL
ED0212XX	Kit, Soft Goods, Athletic, Right Medial / Left Lateral, XS-XXL
ED023YXX	Calf Pad Kit, Everyday
ED025YXX	Calf Pad Kit, Everyday, Short ¹
ED022YXX	Thigh Pad Kit, Everyday
ED024YXX	Thigh Pad Kit, Everyday, Short ¹
ED027YXX	Thigh Pad Kit, Hi-Activity
ED029YXX	Thigh Pad Kit, Hi-Activity, Short ¹

XX = See size chart

Y: Left(1), Right (2)

¹Also available in XXS (XX=00)

Product Sizing Chart

DUO Size Chart				
Size #	Size	Thigh Circumference	Standard Calf Circumference	Athletic Calf Circumference
XX = 01	XS	13.5" - 16" (34 - 41cm)	12.5" - 14" (32 - 36cm)	11" - 12.5" (28 - 32cm)
XX = 03	S	16" - 18.75" (41 - 48cm)	14" - 15.5" (36 - 39cm)	12.5" - 14" (32 - 36cm)
XX = 05	M	18.75" - 21.5" (48 - 55cm)	15.5" - 17" (39 - 43cm)	14" - 15.5" (36 - 39cm)
XX = 07	L	21.5" - 24.25" (55 - 62cm)	17" - 18.5" (43 - 47cm)	15.5" - 17" (39 - 43cm)
XX = 09	XL	24.25" - 27" (62 - 69cm)	18.5" - 20" (47 - 51cm)	17" - 18.5" (43 - 47cm)
XX = 11	XXL	27" - 29.5" (69 - 75 cm)	20" - 21" (51 - 53cm)	18.5" - 20" (47 - 51cm)

Thigh circumference measure 6" (15 cm) above mid-patella. Calf circumference measure 6" (15 cm) below mid-patella.

Resources

- [DUO OA Knee Brace Sell Sheet](#)