



## Roadrunner™ Soft Knee Brace



Sizes: xs, s, m, m+, l, xl, xxl

The Roadrunner provides control and support for patients with ligament injuries and instabilities and is frequently used for post-op treatment and mild OA. Measures 16" in length.

### Common Examples of Use

- Most ACL, PCL, MCL, LCL injuries and instabilities
- Meniscal injuries
- Mild OA

### Product Features

- Available in Airmesh® and neoprene
- Polycentric hinges with adjustable flexion extension stops
- Sleeve and wraparound versions available
- Standard and Post-Op hinge options

Category: Adjustable ROM.

### Instructions for Use

[Instructions for Use](#)

### Part Numbers

## Part Numbers

<b>Part #</b>	<b>Description</b>
1411X	Roadrunner, Airmesh Pull-on
1414X	Roadrunner, Airmesh, Open Back Pull-on
1415X	Roadrunner, Airmesh, Open Back, Wraparound
0694X	Roadrunner, Neoprene Pull-on
0696X	Roadrunner, Neoprene, Open Back, Pull-on
0700X	Roadrunner, Neoprene, Open Back, Wraparound

## Product Sizing Chart

<b>Size #</b>	<b>Size</b>	<b>Thigh Circumference</b>
X=1	XS	12" - 15" (30 - 38 cm)
X=2	S	15" - 18" (38 - 46 cm)
X=3	M	16" - 19" (41 - 48 cm)
X=4	L	18" - 21" (46 - 54 cm)
X=5	XL	20" - 23" (51 - 58 cm)
X=6	XXL	22" - 25" (56 - 64 cm)

Circumference taken 6" (15 cm) above mid patella.

Brace measures 16" (41 cm) in length.