



## Latitude® OA Knee Brace

Sizes: Universal



The Latitude OA knee brace offloads the medial or lateral compartment for patients with mild to moderate unicompartmental osteoarthritis of the knee. The universal design limits inventory needed to cover the patient size range and allows for the brace to be used for both medial and lateral osteoarthritis.

### Common Examples of Use

- Moderate to mild unicompartmental osteoarthritis
- Osteochondral defects
- Articular cartilage defects/repairs

### Product Features

- Single Upright OA brace
- Universal sizing (maximum thigh/calf circumference: 32"/21" or 81.2 cm/53.3 cm)
- Side specific (left or right)
- Offloads both medial and lateral compartments (with the use of a key)
- Lightweight, low profile, and only 14" tall (35.56 cm)
- Easy push button clips with numbered buckles for simple and quick application
- Pivoting straps for patient comfort

Category: Functional OA.

### Part Numbers

<b>Part Number</b>	<b>Item</b>
100630-100	LATITUDE OA KNEE BRACE LEFT
100630-200	LATITUDE OA KNEE BRACE RIGHT
100677-000	LATITUDE OA REPLACEMENT BUCKLES
100678-000	LATITUDE OA REPLACEMENT KEY
100689-100	LATITUDE OA REPLACEMENT PAD KIT LT
100689-200	LATITUDE OA REPLACEMENT PAD KIT RT
100690-000	LATITUDE OA REPLACEMENT STRAP KIT
100691-000	LATITUDE OA SUSPENSION SLEEVE

## **Instructions for Use**

[Instructions for Use](#)